

A close-up photograph of a dental light fixture. The fixture is white and circular, with two small rectangular slots on its face. It is mounted on a blue, curved metal arm. The background is dark blue, and the lighting is soft, highlighting the texture of the fixture and the metallic sheen of the arm.

Dental Light
701

Designing a new dental light with good optical capability, maneuverability and design with reasonable cost is the basic concept behind 701



Easy to replace new halogen lamp

The three-position controller enables the selection of composite safety mode to prevent premature curing and two intensity modes, low and high for proper oral cavity illumination

Comfort grip handles facilitate preferred positioning

Compact and Sleek design light head

Touchless sensor activates the hands free on and off operation, offering reduced risk of cross contamination

Options

Detachable Handles



Autoclavable and easy to remove.

Mirror



Adjustable angle of the patient mirror to look inside mouth.

Manual Switch



Manual Switches for on-off and three-positions control, low/high and composites mode.



**Unit Mount Light
Model; 720S**

The unit mount lights offer a great option for smaller surgeries.

**Ceiling Mount Light
Model; 702S**

The ceiling mount lights offer a wider space for operating area.



**Track Light
Model; 7055**

The track lights offer additional flexibility for larger surgeries.

