


Dental Chair

PRO |

037/039



 **Belmont**[®]

Protection for Professionals

PRO ■ O37/O39

O39

AXIS POWER HEADREST

moves smoothly due to its axis mechanism which simulates rotation of patient's neck, and patient can comfortably keep mouth open and provides dentist with an improved view of oral cavity.




O37

OFFSET ARTICULATING HEADREST

provides occipital support and proper angulation for viewing oral cavity.

(OPTIONAL FEATURE FOR O37) TWIN AXIS ARTICULATING HEADREST

is reversible and suitable for both children and adults.



SWING-OUT ARMREST

for easy entry and exit for patient.

LUMBER SUPPORT

is anatomically correct in all positions.

UPHOLSTERY

Easy to clean, seamless upholstery.

Flame retardant.

A choice from 11 colours.



PRO I CHAIR FEATURES


This unique design create; communication between <

HANDS-FREE OPERATION

3 preset positions, back £
Dentist, staff, and patient &
contamination.

Eye to oral cavity contact

THE SYNCHRONOUS MO



ROTATION

Chair rotates to maximum
of 110° to Left & Right.

FOOT CONTROLS

allow operator to position patient quickly and
correctly leaving hands free and clean.



A FOLDING LEGREST WHICH PROVIDES EASY ENTRY AND EXIT
an ideal environment for personal consultation and
dentist and patient.

PERMITS EFFECTIVE INFECTION CONTROL.

and base movements are operated by foot lever.

•e protected from disease transmission from finger-to-switch

can be maintained during procedure.

EMENT OF BACKREST AND SEAT PERFECTLY FITS PATIENT'S BODY


STATUS INDICATOR PANEL

incorporates simple to operate chair position
programmes (max.8) and self diagnostic function.



ROLL-UP LEGREST

extends when chair is in a preset position
and retracts when chair is in a consultation
position or entry/exit position.



O39

FOOT SWITCH FOR HEADREST

Hands-free operation of power headrest
eliminates risk of cross contamination
from hand-to-headrest contact.




- 1) LOCK PEDAL
- 2) MANUAL (BACK & BASE MOVEMENTS)
- 3) AUTO-RETURN/2 PRESETS/LAST POSITION MEMORY
- 4) MANUAL (BACK & BASE MOVEMENTS]

PRO



Dimensions/Specifications



The manufacturers reserve the right to alter and improve the specification without notice.

Printed in Japan 2000.05MA